

cabinet-api-0.1

CabiNET API

API Version: 0.1

INDEX

1. API	4
1.1 GET /api/v1/device/list	4
1.2 POST /api/v1/device/provision/start	4
1.3 POST /api/v1/device/provision/{id}/verify	5
1.4 GET /api/v1/device/{id}	5
1.5 PUT /api/v1/device/{id}	6
1.6 GET /api/v1/device/{id}/dispense_time	6
1.7 POST /api/v1/device/{id}/dispense_time	7
1.8 GET /api/v1/device/{id}/medication	8
1.9 POST /api/v1/device/{id}/medication	8
1.10 GET /api/v1/device/{id}/medication/{med_id}	9
1.11 DELETE /api/v1/device/{id}/medication/{med_id}	10
1.12 POST /api/v1/device/{id}/reload	10
1.13 POST /api/v1/device/{id}/state	11
1.14 POST /api/v1/device/{id}/sync	11
1.15 POST /api/v1/user/anonymous_upgrade	12
1.16 GET /api/v1/user/me	12
1.17 POST /api/v1/user/register	12
1.18 POST /api/v1/user/register_anonymous	13
1.19 POST /api/v1_2/device/auth	13
1.20 POST /api/v1_2/device/provision	13
1.21 POST /api/v1_2/device/sync	14
2. JWK-SET.JSON	15
2.1 GET /jwk-set.json	15
3. SECURITY	16
3.1 POST /api/v1/auth/login	16
3.2 POST /api/v1/auth/login_anonymous	16

Security and Authentication

SECURITY SCHEMES

KEY	TYPE	DESCRIPTION
jwt	http, bearer, jwt	

API

1. API

1.1 GET /api/v1/device/list

Lists devices that the user has access to

REQUEST

No request parameters

RESPONSE

STATUS CODE - 200: listDeviceUser 200 response

RESPONSE MODEL - application/json

```
[{  
  Array of object:  
    id*           integer  
    deviceID*     integer  
    deviceClass*  enum      ALLOWED:v1_7x2  
                                Describes the "type" of pill organizer a device is. There is presently only a single pill organizer  
                                model. This enum is to possibly support multiple models in the future, possibly with wildly  
                                different bin configurations.  
    customName    string  
    lastSync      string  
    serialNo*     integer  
    primaryUser*  boolean  
    owner*        boolean  
    notifications* boolean  
    timezone      string  
}]
```

1.2 POST /api/v1/device/provision/start

Initiates provisioning

REQUEST

REQUEST BODY - application/json

```
{  
  serialNo  integer  
  deviceClass string  
}
```

HEADER PARAMETERS

NAME	TYPE	DESCRIPTION
X-Local-TZ	string	

RESPONSE

STATUS CODE - 200: startProvisioning 200 response

RESPONSE MODEL - application/json

```
{
  A record of a particular attempt at provisioning a pill organizer device, regardless of it was successful or not. This is essentially a "session"
  record of a provision attempt. A device can potentially have many different provision records.

  id          integer
  active       boolean
  oobKey       [string]
  userID       integer
  timezone     string
}
```

1.3 POST /api/v1/device/provision/{id}/verify

Checks provisioning status

REQUEST

QUERY PARAMETERS

NAME	TYPE	DESCRIPTION
*id	int64	

REQUEST BODY - application/json

```
{
  serialNo integer
  ssid     string
}
```

RESPONSE

STATUS CODE - 200: checkProvisionStatus 200 response

RESPONSE MODEL - application/json

```
{
}
```

1.4 GET /api/v1/device/{id}

Queries info about a device

REQUEST

QUERY PARAMETERS

NAME	TYPE	DESCRIPTION
*id	int64	

RESPONSE

STATUS CODE - 200: device 200 response

RESPONSE MODEL - application/json

```
{
  A pill organizer device.

  id          integer
  userType     enum      ALLOWED:STANDARD, DEVICE, ANONYMOUS
}
```

```
deviceClass  enum
    serialNo   integer
    customName string
    lastSync   string
    ipv4        integer
    ipv6        [string]
    baseTZ      string
}
```

The purpose of this enumeration was to differentiate between devices (which can login) and regular users. This is a bad idea and must be refactored out.

ALLOWED:v1_7x2

Describes the "type" of pill organizer a device is. There is presently only a single pill organizer model. This enum is to possibly support multiple models in the future, possibly with wildly different bin configurations.

1.5 PUT /api/v1/device/{id}

Updates device basic settings

Updates non-schedule settings for a device, including timezone, name, and notification token.

REQUEST

PATH PARAMETERS

NAME	TYPE	DESCRIPTION
*id	int64	

REQUEST BODY - application/json

```
{
  deviceName*    string
  notificationToken* string
  notifications*  boolean
  timezone*      string
}
```

RESPONSE

STATUS CODE - 200: setDeviceSettings 200 response

RESPONSE MODEL - application/json

```
{
  id*          integer
  deviceID*    integer
  deviceClass* enum    ALLOWED:v1_7x2
                  Describes the "type" of pill organizer a device is. There is presently only a single pill organizer model.
                  This enum is to possibly support multiple models in the future, possibly with wildly different bin
                  configurations.

  customName    string
  lastSync      string
  serialNo*     integer
  primaryUser*  boolean
  owner*        boolean
  notifications* boolean
  timezone      string
}
```

1.6 GET /api/v1/device/{id}/dispense_time

Get dispense times

REQUEST

PATH PARAMETERS

NAME	TYPE	DESCRIPTION
*id	int64	

RESPONSE

STATUS CODE - 200: dispenseTimes 200 response

RESPONSE MODEL - application/json

```
{
  A simplified representation of a "DeviceSimpleScheduleStrategy" and two "DeviceSimpleDispenseTime" objects, hiding the underlying complexity of the data structure.
  amID integer ID of the AM dispense time
  amSecondsFrom00 integer AM dispense time in day-epoch format (seconds)
  pmID integer ID of the PM dispense time
  pmSecondsFrom00 integer PM dispense time in day-epoch format (seconds)
}
```

1.7 POST /api/v1/device/{id}/dispense_time

Update dispense times

REQUEST

PATH PARAMETERS

NAME	TYPE	DESCRIPTION
*id	int64	

REQUEST BODY - application/json

```
{
  A simplified representation of a "DeviceSimpleScheduleStrategy" and two "DeviceSimpleDispenseTime" objects, hiding the underlying complexity of the data structure.
  amID integer ID of the AM dispense time
  amSecondsFrom00 integer AM dispense time in day-epoch format (seconds)
  pmID integer ID of the PM dispense time
  pmSecondsFrom00 integer PM dispense time in day-epoch format (seconds)
}
```

RESPONSE

STATUS CODE - 200: updateDispenseTime 200 response

RESPONSE MODEL - application/json

```
{
  A simplified representation of a "DeviceSimpleScheduleStrategy" and two "DeviceSimpleDispenseTime" objects, hiding the underlying complexity of the data structure.
  amID integer ID of the AM dispense time
  amSecondsFrom00 integer AM dispense time in day-epoch format (seconds)
  pmID integer ID of the PM dispense time
  pmSecondsFrom00 integer PM dispense time in day-epoch format (seconds)
}
```

1.8 GET /api/v1/device/{id}/medication

List medications for a device

REQUEST

PATH PARAMETERS

NAME	TYPE	DESCRIPTION
*id	int64	

RESPONSE

STATUS CODE - 200: medications 200 response

RESPONSE MODEL - application/json

```
[{
  id integer
  dispenseTimes [{
    id integer
    dispense {
      id integer
    }
    medicationID integer
    dispenseID integer
    quantity integer
  }]
  deviceID integer
  med_name* string 1 to 64 chars
  shape* enum ALLOWED:CAPSULE, DIAMOND, DOUBLE_CIRCLE, HEPTAGON, HEXAGON,
  OCTAGON, OVAL, PENTAGON, RECTANGLE, ROUND, SEMICIRCLE, SQUARE,
  TEAR, TRAPEZOID, TRIANGLE
  color integer
}]
```

1.9 POST /api/v1/device/{id}/medication

Saves a medication

REQUEST

PATH PARAMETERS

NAME	TYPE	DESCRIPTION
*id	int64	

REQUEST BODY - application/json

```
{
```



```

id                integer
name*             string
shape*            enum      ALLOWED:CAPSULE, DIAMOND, DOUBLE_CIRCLE, HEPTAGON, HEXAGON,
                             OCTAGON, OVAL, PENTAGON, RECTANGLE, ROUND, SEMICIRCLE,
                             SQUARE, TEAR, TRAPEZOID, TRIANGLE

color*            integer
dispenseTimes*    [integer]
}

```

RESPONSE

STATUS CODE - 200: saveMedication 200 response

RESPONSE MODEL - application/json

```

{
  A medication that is taken according to some schedule. We currently view a "medication" to be anything you could put in a pill organizer - currently
  we only really consider pills under this view. A `ScheduledMedication` has a name, shape, color, and a frequency and periodicity.

  id                integer
  dispenseTimes [ {
    Array of object: Relates a `ScheduledMedication` to a `DeviceBaseDispenseTime`. Also contains the quantity of medication to be taken
    at a given dispense time.

    id                integer
    dispense {
      A single medication dispense time, for a particular device's schedule. Medications are attached to a dispense time to indicate they
      should be taken then.

      id integer
    }
    medicationID      integer
    dispenseID        integer
    quantity           integer
  } ]
  deviceID            integer
  med_name*           string  1 to 64 chars
  shape*              enum    ALLOWED:CAPSULE, DIAMOND, DOUBLE_CIRCLE, HEPTAGON, HEXAGON,
                              OCTAGON, OVAL, PENTAGON, RECTANGLE, ROUND, SEMICIRCLE, SQUARE,
                              TEAR, TRAPEZOID, TRIANGLE

  color              integer
}

```

1.10 GET /api/v1/device/{id}/medication/{med_id}

Gets a medication by ID

REQUEST

PATH PARAMETERS

NAME	TYPE	DESCRIPTION
*id	int64	
*med_id	int64	

RESPONSE

STATUS CODE - 200: get 200 response

RESPONSE MODEL - application/json

```

{
  A medication that is taken according to some schedule. We currently view a "medication" to be anything you could put in a pill organizer - currently
  we only really consider pills under this view. A `ScheduledMedication` has a name, shape, color, and a frequency and periodicity.

  id integer
  dispenseTimes [{
    Array of object: Relates a `ScheduledMedication` to a `DeviceBaseDispenseTime`. Also contains the quantity of medication to be taken
    at a given dispense time.

    id integer
    dispense {
      A single medication dispense time, for a particular device's schedule. Medications are attached to a dispense time to indicate they
      should be taken then.

      id integer
    }
    medicationID integer
    dispenseID integer
    quantity integer
  }]
  deviceID integer
  med_name* string 1 to 64 chars
  shape* enum ALLOWED:CAPSULE, DIAMOND, DOUBLE_CIRCLE, HEPTAGON, HEXAGON,
  OCTAGON, OVAL, PENTAGON, RECTANGLE, ROUND, SEMICIRCLE, SQUARE,
  TEAR, TRAPEZOID, TRIANGLE
  color integer
}

```

1.11 DELETE /api/v1/device/{id}/medication/{med_id}

Deletes a medication

REQUEST

PATH PARAMETERS

NAME	TYPE	DESCRIPTION
*id	int64	
*med_id	int64	

RESPONSE

STATUS CODE - 200: deleteMedication 200 response

1.12 POST /api/v1/device/{id}/reload

Reloads the pills on a device

REQUEST

QUERY PARAMETERS

NAME	TYPE	DESCRIPTION
*id	int64	

RESPONSE

STATUS CODE - 200: reload 200 response

1.13 POST /api/v1/device/{id}/state

Get device state on a particular date

REQUEST

PATH PARAMETERS

NAME	TYPE	DESCRIPTION
*id	int64	

QUERY PARAMETERS

NAME	TYPE	DESCRIPTION
*date	string	

RESPONSE

STATUS CODE - 200: consolidatedState 200 response

RESPONSE MODEL - application/json

```
{
  id*           integer
  lastSync      string
  bins*         integer
  dosePeriods* [{
    Array of object:
      binID*     integer
      timestamp* integer
      status*    integer
      medications* [integer]
  }]
}
```

1.14 POST /api/v1/device/{id}/sync

App-proxied device sync

REQUEST

PATH PARAMETERS

NAME	TYPE	DESCRIPTION
*id	int64	

REQUEST BODY - application/json

string

RESPONSE

STATUS CODE - 200: sync 200 response

RESPONSE MODEL - application/json

```
{  
}
```

1.15 POST /api/v1/user/anonymous_upgrade

Upgrades an anonymous account into a standard account

Performs an in-place upgrade of an anonymous account into a standard account with an email and password. The user ID is preserved.

REQUEST

REQUEST BODY - application/json

```
{  
  email    string  
  password string 6 to 32 chars  
}
```

RESPONSE

STATUS CODE - 200: upgradeAnonymous 200 response

RESPONSE MODEL - application/json

```
{  
  A standard user account with an email/password combo. Users of this type are intended to be humans who explicitly registered for a CabiNET  
  account and voluntarily supplied their own email and password.  
  id      integer  
  role    enum      ALLOWED:USER, ADMIN  
           An enumeration of possible special roles for a user.  
  email   string  
}
```

1.16 GET /api/v1/user/me

Gets info about currently signed-in user

REQUEST

No request parameters

RESPONSE

STATUS CODE - 200: authenticationStatus 200 response

RESPONSE MODEL - application/json

```
{  
  id*    integer  
  email* string  
}
```

1.17 POST /api/v1/user/register

Register a standard account.

REQUEST

REQUEST BODY - application/json

```
{
```

```
email    string
password string 6 to 32 chars
}
```

RESPONSE

STATUS CODE - 200: register 200 response

RESPONSE MODEL - application/json

```
{
  A standard user account with an email/password combo. Users of this type are intended to be humans who explicitly registered for a CabiNET
  account and voluntarily supplied their own email and password.
  id      integer
  role    enum      ALLOWED:USER, ADMIN
              An enumeration of possible special roles for a user.
  email   string
}
```

1.18 POST /api/v1/user/register_anonymous

Creates an anonymous user

REQUEST

No request parameters

RESPONSE

STATUS CODE - 200: anonymousRegistration 200 response

RESPONSE MODEL - application/json

```
{
  An `AnonymousUser` is a user who has not formally registered for an account but still needs to be associated with some sort of user for
  authentication and authorization purposes. This allows customers to make use of the app without having to sign up for a full account. The
  `secret` is designed to be stored in secure on-device storage and not shared with the user. An anonymous user can be upgraded to a full user.
  id      integer
  role    enum      ALLOWED:USER, ADMIN
              An enumeration of possible special roles for a user.
  secret  [string]
}
```

1.19 POST /api/v1_2/device/auth

Authenticates a device using their credentials in a Protobuf AuthorizeRequest

REQUEST

RESPONSE

STATUS CODE - 200: login_1 200 response

RESPONSE MODEL - application/x-protobuf

```
{
}
```

1.20 POST /api/v1_2/device/provision

Completes the provisioning process

Indicates that the device has been successfully provisioned and connected to WiFi. The device should update its state using the sync response provided by this endpoint.

REQUEST

RESPONSE

STATUS CODE - 200: provision 200 response

RESPONSE MODEL - application/x-protobuf

```
{  
}
```

1.21 POST /api/v1_2/device/sync

Syncs a device's state

Performs a two-way device sync, accepting a device's state, processing events, and returning the server's authoritative state.

REQUEST

RESPONSE

STATUS CODE - 200: sync_1 200 response

RESPONSE MODEL - application/x-protobuf

```
{  
}
```

2. JWK-SET.JSON

2.1 GET /jwk-set.json

REQUEST

No request parameters

RESPONSE

STATUS CODE - 200: getJWKS 200 response

RESPONSE MODEL - application/json

```
{  
}
```

3. SECURITY

3.1 POST /api/v1/auth/login

accepts an email/password pair to authenticate a JWT session

REQUEST

REQUEST BODY - application/json

```
{
  username string The account's email address
  password string Password to authenticate the account with
}
```

FORM DATA PARAMETERS

NAME	TYPE	DESCRIPTION
username	string	
password	string	
identity	string	
secret	string	

RESPONSE

STATUS CODE - 200: authentication success

RESPONSE MODEL - application/json

```
{
  username      string  The account's email address
  access_token  string  JWT token
  token_type    string  The type of JWT token (this API supports only Bearer tokens)
  expires_in    integer The amount of seconds the newly-issued JWT token is valid for
}
```

STATUS CODE - 401: authentication failure

RESPONSE MODEL - application/json

```
{
  message string A short message describing why authentication failed
}
```

3.2 POST /api/v1/auth/login_anonymous

accepts an id/secret pair

REQUEST

REQUEST BODY - application/json

```
{
  id      number ID
  secret  string Secret
}
```

RESPONSE

STATUS CODE - 200: authentication success

RESPONSE MODEL - application/json

```
{
  username      string  The account's email address
  access_token  string  JWT token
  token_type    string  The type of JWT token (this API supports only Bearer tokens)
  expires_in    integer The amount of seconds the newly-issued JWT token is valid for
}
```

STATUS CODE - 401: authentication failure

RESPONSE MODEL - application/json

```
{
  message string A short message describing why authentication failed
}
```
